RVENT

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BREATHING SYSTEMS INSTRUCTION FOR USE

The medical device meets the requirements of the ISO 80601-2-13 and ISO 5356-1 standards.

INTENDED USE:

Anesthesia and Ventilation Circuits: Disposable breathing circuit for conduction of respiratory gases between anesthesia machine or ventilator and patient and intended for single use only. Sterile and Non-sterile options are available.

Breathing bag with connection hose (limb) intended for use with anesthesia delivery systems as a reservoir during automatic ventilation and as a manual breathing bag during manual ventilation.

Heated Wire Breathing Circuits: Disposable heated wire breathing circuit for conveying moistened breathing gas between the humidifier and patients. Intended for single use only. Sterile and Non-sterile options are available.

IPPB&Homecare Breathing Circuits: The IPPB breathing systems are designed to be used with Intermittent Positive Pressure Breath machines for spontaneously breathing patients, to assist lung expansion, deliver aerosol therapy or to assist ventilation. Sterile and Non-sterile options are available.

BPAP Circuits: Bilevel Positive Airway Pressure (BPAP) is a modified version of the continuous positive airway pressure (CPAP). It is a mode of ventilation with cycling variations between two continuous positive airway pressure levels. This mode offers two different levels of air pressure: Increased pressure during inhalation and decreased pressure during exhalation. Sterile and Non-sterile options are available.

CPAP Circuits: Continuous Positive Airway Pressure (CPAP) is the method of ventilating a patient applying higher pressure than environmental atmospheric pressure to the airways of the patient to improve oxygenation via ventilator or CPAP machine. Sterile and Non-sterile options are available.

WARNINGS:

General Warnings:

- The medical device must only be used as specified in the intended use section.
- The medical device manufactured for single use only.
- Use the breathing circuit only with components contained in the set.
- The medical device must not be reused or reprocessed. Reusing or reprocessing of the product may cause the patient injury. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.
- To avoid contamination and soiling, keep the medical device packaged until ready to be used. Do not use the medical device if the packaging is damaged.
- Do not modify the medical device. Modification may damage or impair the proper functioning of the device which may lead to patient injury.
- Obstructions, damage, and foreign matter can lead to malfunction. Check all system components for obstructions, damage, and foreign matter before installation.
- Check the breathing set for damage prior to use. If there is any damage, do not use the respirator and contact the manufacturer.
- Ensure that all connections are secure and free of leakages.
- Perform a selftest of the basic device including a leakage test after the breathing circuit (hose, filter/HME, etc.) has been completely installed and before use on the patient.
- During installation position the breathing circuit so that the risk of tripping or of stepping on the breathing circuit is minimized.
- Make sure to install the breathing circuit without loops and kinks as they may increase resistance.
- If too much condensate accumulates, a partial or complete blockage of the breathing circuit can occur. Regularly check the breathing circuit for condensate and empty if necessary.
- Patient parameters, the patient's clinical status, and the properties of components of the breathing circuit (e.g., tube) may influence or limit the performance characteristics of the breathing circuit during high-frequency ventilation. Ensure appropriate patient monitoring.
- Adding further components and using incompatible components may
- increase inspiratory and expiratory resistance and adversely affect the performance of the ventilator.
- If includes a water trap; wrong installation can lead to malfunction of the
- water trap. Install the water trap correctly. If the water trap is installed higher than the patient, condensate can flow into the breathing hoses. Install the water trap at the lowest point of the breathing circuit and below patient level.
- If includes; A Luer Lock connector must only be used for gas monitoring. Any other use of the Luer Lock connector may endanger the patient. When not in use, the Luer Lock connector must be sealed with its cap.
- If the circuit is coaxial, before using the device on the patient, check the inspiratory line for leaks in the inner tubing using a leakage test device.

If the product is packaged sterile;

- The device is sterilized with Ethylene Oxide. Sterile conditions are guaranteed only if the package is not damaged.

STERILE EO

- Do not re-sterilize. Re-sterilizing damages the product integrity.

Additional Warnings for Heated Wire Breathing Circuits:

- The heated breathing circuit must never be operated without gas flow. The humidifier must be switched off if the gas flow is interrupted.
- Fill the humidifier chamber with distilled water before using the breathing circuit and check the function of the floater.
- Only use distilled water. The use of other substances can harm the patient.
- The humidifier chamber must be positioned lower than the patient.
- When using active humidification, do not use the breathing circuit in conjunction with a heat and moisture exchanger (HME). There is a risk of pressure buildup and of insufficient ventilation as a result of water accumulating in the HME.
- To prevent from risk of burns, avoid prolonged and direct contact between the breathing circuit and the patient's skin.
- To prevent from risk of burns, the heated breathing circuit must not be covered with bed sheets, hand towels or similar flammable objects.
- Too high or too low ambient temperatures can lead to increased accumulation of humidity or humidification performance can be reduced.
- Do not operate the device if the water level in the humidifier chamber exceeds the marked maximum.
- Do not touch the hot plate or the bottom of the chamber. The surfaces can reach a temperature of over 85°C. These hot surfaces radiate heat.

Protection Class of Heated Wire Breathing Circuits: Type BF

INSTRUCTION FOR USE:

Anesthesia and Ventilation Circuits: Connect the inspiration input connector to the appropriate inspiration port specified on the device. Connect the expiration input connector to the expiration port specified in the device. If the circuit contains a limb, connect the limb to the port on the device in a suitable way. If the circuit contains breathing bag, place the breathing bag on the end of the limb. If the circuit contains a gas sampling line, connect one end of the gas sampling line to the port on the elbow connector and the other to the capnograph device. If the circuit contains a breathing filter, place it in the appropriate position. If the circuit contains a water trap, make sure that is fully closed. During the draining of the water contained in the water trap, turn it counterclockwise and open it to drain it. If the circuit contains an breathing mask, adjust it through the valve on the cushion with the syringe to bring the cushion to the desired level and connect it with the circuit.

Heated Wire Breathing Circuits: Insert the humidifier chamber completely into the humidifier until it clicks into place. To remove it, pull the humidifier chamber out of the humidifier. Connect the ventilator inspiratory limb first to the humidifier, and then to the inspiratory port on the ventilator. Connect the ventilator expiratory limb to the expiratory port on the ventilator and the Y-piece. Connect the patient inspiratory limb first to the humidifier chamber, and then to the Y-piece. Insert the humidifier chamber water feed spike into the water reservoir. Note that the limb connectors on the humidifier combine electrical connections with breathing circuit connectors. Ensure proper orientation of electrical contacts on the breathing circuit connectors to match the connecting element on the humidifier.

IPPB&Homecare Breathing Circuits: Unpack the product. Attach the circuit to IPPB transport ventilator. Pressure test the circuit and machine to insure proper function. Instruct patient to purse lips around mouthpiece (if includes) so air do not leak. Monitor patient throughout duration of treatment. After the treatment is completed detach IPPB circuit from the patient.

BPAP/CPAP Breathing Circuits: Unpack the product. Attach the one end of the circuit to ventilator or CPAP/BPAP machine. Attach the other end of the circuit to the mask. Ensure proper fitting. After the treatment is completed, detach the circuit from the patient.

INDICATIONS:

Anesthesia and Ventilation Breathing Circuits: The circuits are designed to connect the patient to the respiratory gas supply during anesthesia and intensive care.

Heated Wire Breathing Circuits: It provides heating and delivery of respiratory gases to patients who need respiratory support.

IPPB&Homecare Breathing Circuits: The IPPB breathing circuits promote the facilitation of excretion of secretions, the expansion of the bronchial lumen, increased ventilation of the alveoli and better oxygen passage into the blood in chronic care patients.

BPAP Breathing Circuits: BPAP is used for patients with restrictive or obstructive disorders and for patients with sleep apnea who can not tolerate the high pressure of CPAP.

CPAP Breathing Circuits: Obstructive sleep apnea syndrome, Acute respiratory distress syndrome, Post-operative atelectasis, Acute respiratory distress, Chronic obstructive pulmonary disease, Cardiogenic pulmonary oedema, Post-extubation hypoxemia, Weaning.

CONTRAINDICATIONS:

Anesthesia and Ventilation Breathing Circuits: The device is used only as directed. Improper use of it will harm the patient.

Heated Wire Breathing Circuits: The device must be used only as directed. Improper use of it will harm the patient.

IPPB&Homecare Breathing Circuits: Increased intracranial pressure, Hiccups, Hemodynamic instability, Recent facial, oral or skull surgery, Tracheoesophageal fistula, Recent esophageal surgery, Active haemoptysis, Nausea, Active untreated tuberculosis or other respiratory communicable disease, Radiographic evidence of bleb.

BPAP/CPAP Breathing Circuits: The device must be used only as directed. Improper use of it will harm the patient.

TIME OF USAGE:

All of the Breathing Systems are intended for single use only. It is used during the application period.

SHELF LIFE:

For all of the Breathing Systems, Shelf Life is 5 years.

TARGET PATIENT GROUP:

Adult, pediatric and neonatal patients.

USERS:

Users have the necessary expertise, training, and knowledge to use or install the medical device.

RISK OF FIRE:

In combination with oxygen or nitrous oxide, ignition sources such as electrosurgery and laser surgery devices can cause fires. To protect patient and users, prevent leakages from hose carrying oxygen or nitrous oxide.

In order not to damage the breathing circuit, keep a distance of at least 200 mm (7.9 in) between hoses carrying oxygen or nitrous oxide and a possible ignition source (e.g., electrosurgery or laser surgery devices).

DISPOSAL METHOD:

Local regulations and/or hospital waste management procedures of the relevant country should be followed when disposing of the used products in order to avoid possible contamination.

PERFORMANCE DATA:

The performance data conforms to the relevant sections of the ISO 5367 standard. All breathing tube types (corrugated, smoothbore and extendible) and patients groups (adult, pediatric and neonatal) tested for resistance to flow and compliance according to relevant sections of ISO 5367 standard.

Resistance to Flow at 30L/min:	<0.06hPa//min
Resistance to Flow at 15L/min:	<0.12hPa//min
Resistance to Flow at 2.5L/min:	<0.74hPa//min
Leakage at 60mbar:	For Adult: <70ml/min
	For Pediatric: <40ml/min
	For Neonatal: <30ml/min
Compliance at 60 hPa:	For Adult: <5ml/hPa
	For Pediatric: <4ml/hPa
	For Neonatal: <1.5ml/hPa

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